

CONTACT US NOW

+44 (0) 1432 850 910

support@medDXsolutions.com

www.medDXsolutions.com

# SAP Material Super Absorbent Pads


## Products

## Components

## Sample Collection Kits

## Guidance

Product Code	Case Qty.	Description (blood absorbent capacity)	Dimensions
AB001/1000	1000	SAP 50ml x 1000	75 x 125 mm
AB001/500	500	SAP 50ml x 500	75 x 125 mm
AB001/100	100	SAP 50ml x 100	75 x 125 mm
AB002/500	500	SAP 100ml x 500	125 x 150 mm
AB002/100	100	SAP 100ml x 100	125 x 150 mm
AB003/500	500	SAP 200ml x 500	150 x 250 mm
AB003/100	100	SAP 200ml x 100	150 x 250 mm
AB004/250	250	SAP 400ml x 250	250 x 300mm
AB004/100	100	SAP 400ml x 100	250 x 300mm
AB010/500	500	SAP 4 Sleeve x 500	
AB010/100	100	SAP 4 Sleeve x 100	



Produced from a Super Absorbent Polymer material, MedDX SAP sample packing absorbent sheets are designed to help meet UN3373 PI650 regulations by providing excellent absorption capacity. The size of sheet used should be that whose absorbent capacity at least matches the total volume of liquid being shipped in the primary containers (see table).

MedDX SAP absorbent is also available as a 4 Bay sleeve, allowing for multiple samples to be separated and protected.

**MedDX**Solutions  
A More Complete Approach